

RenMab

An Innovative Platform of Biocytogen




Key Features:

- Developed using Biocytogen's unique Mb-scale chromosome engineering technology
- Variable regions of mouse heavy chain and κ light chain genes are replaced *in situ* by human counterparts
- Constant region is retained to ensure proper B cell development
- Normal immune organs development and good immune response to immunization like WT mouse
- RenMab antibodies exhibit human-like repertoire and high affinity to antigens in subnanomolar range

Generation of RenMab™ Mouse


Biocytogen has developed a fully human antibody mouse (RenMab™ Mouse), whose genes that encode entire antibody variable regions are replaced *in situ* by human Ig heavy chain and κ light chain through Mb-scale chromosome engineering technology. RenMab™ Mouse provides an efficient platform for fully human antibody generation, characterization, therapeutic antibody discovery, and rapid *in vivo* efficacy screening.

Mouse variable regions of the heavy and κ light chains are replaced by the human counterparts *in situ*



RenMab™ Mouse demonstrates a similar immune system to that of WT mouse

RenMab Mice immune cells profiling suggest a comparable immune system with wild type mice




Full human VDJ expression in the variable domains of heavy chain and kappa light chain genes


The RenMab™ mouse shows WT-like expression of antibody levels and full usage of human heavy chain V, D, and J domain segments. Similarly, mouse kappa chain segments were replaced with human V and J segments and also show WT-like expression of antibody levels.

RenMab™ Naive Mouse Heavy Chain IGHV Germline Usage


RenMab™ Naive Mouse Heavy Chain IGHD & IGHJ Germline



RenMab™ Naive Mouse Kappa Chain IGKV Germline Usage




RenMab™ Naive Mouse Kappa Light Chain IGKJ Germline Usage




RenMab™ Mouse produces human-like antibody repertoire as shown by CDR3 analysis

CDRH3 length (amino acids) IgM-Native RenMab™ spleen cell (n=2)



Heavy chain CDR3 length distribution: median length 14.2 amino acids


RenMab™ Heavy Chain CDR3 Amino Acid Usage Frequency



Human-like heavy chain V-D-J recombination pattern: (1) length-dependent usage of tyrosine; (2) preferred usage of D, E and G at IMGT position 107; (3) IMGT position 115 and 117 are conserved due to human JH4 and JH6 usage; (4) Usage of cysteine residue in CDR3, indicative of human DH2 gene family usage in V-D-J recombination


RenMab™ Mouse demonstrates highly similar immune responses compared to WT mouse

Robust immune response in RenMab™ elicited by a panel of antigens



- Robust immune responses against strong/weak antigens were observed
- Accelerated immunization strategies are being evaluated
- Several immunization strategies for difficult targets are being optimized for better outcome
- RenMab-based KO antibody discovery is ongoing


Consistent Ig isotype ratios in naive and immunized RenMab™ mice suggested successful class switch of antibodies



No significant differences in serum total IgG1, IgG2b, IgG2c, IgG3 and IgM levels were observed between RenMab and C57BL/6

RenMab™ antibodies exhibit high affinity against antigens at subnanomolar range

Clone ID	ka (1/Ms)	kd (1/s)	KD (M)	Rmax (RU)
03-7C7	4.35E+05	1.06E-04	2.43E-10	60.1
03-10D12	3.56E+05	1.07E-04	3.00E-10	75.8
09-4E5	3.54E+05	1.17E-04	3.30E-10	98.3
03-9A8	2.35E+05	9.86E-05	4.20E-10	68.0
03-3E7	4.72E+05	2.46E-04	5.22E-10	61.8
02-5A8	3.20E+05	1.68E-04	5.24E-10	56.5
05-8A10	3.19E+05	2.38E-04	7.44E-10	68.7
03-6E12	4.47E+05	3.63E-04	8.14E-10	58.5
05-8G6	3.37E+05	2.81E-04	8.32E-10	85.9
03-1F9	2.00E+05	1.67E-04	8.35E-10	90.9
02-5C1	3.67E+05	3.63E-04	9.90E-10	34.6
02-4B5	4.16E+05	4.25E-04	1.02E-09	37.6
03-9H12	4.26E+05	4.98E-04	1.17E-09	54.5
05-9G11	5.89E+05	6.87E-04	1.17E-09	62.5
09-7B8	1.86E+05	3.13E-04	1.68E-09	24.7
03-10C10	2.74E+05	4.74E-04	1.73E-09	4.9
05-10B7	4.99E+05	9.09E-04	1.82E-09	38.6
02-5C1	4.34E+05	9.54E-04	2.20E-09	31.9
02-5A8	3.27E+05	7.59E-04	2.32E-09	30.3
01-7G7	3.80E+05	9.23E-04	2.43E-09	68.3
05-7C1	1.91E+05	5.75E-04	3.01E-09	55.6



RenMab Commercial Benefits




Licensing Options

License directly with Biocytogen
Option to use through other CROs or in-house



Co-development Opportunities

Exclusive partnership and co-development opportunities for Project Intergrum



Flexibility

Flexible terms tailored to accommodate different antibody programs

CONTACT US



50C Audubon Road
Wakefield, MA 01880, USA



781-587-3558



info@biocytogen.com



www.renmab.com



BIOCYTOGEN